



JHARKHAND
Rai University

Established by an Act of Govt. of Jharkhand
as per Section 2f of UGC Act, 1956

PRACTICAL LAB MANUAL

LAB MANUAL PHARMACEUTICAL

ANALYSIS - i

(B. Pharm Ist Year)

TABLE OF CONTENT

S. No.	Name of Experiment	Date	Remarks	Signature of Teacher
1.	To study the steps involved in volumetric analysis.			
2.	To prepare and standardize 0.1 N sodium hydroxide solution (100 ml).			
3.	To perform the assay of ammonium chloride by acid base titration.			
4.	To prepare & standardize 0.1 M sodium thiosulphate solution (100 ml).			
5.	To prepare & standardize 0.02 N potassium permanganate solution (100 ml).			
6.	To perform the assay of the dilute hydrogen peroxide solution by permanganometry.			
7.	To determine the normality of strong acid by conductometric titration.			
8.	To determine the normality of strong base by potentiometric titration.			
9.	To prepare and standardize 0.05 M disodium edetate solution (100 ml).			
10.	To draw neutralization curve between the given acid-base samples & predict the indicator\s that can be used.			

EXPERIMENT NO. 1

1. OBJECTIVE

To study the steps involved in volumetric analysis.

2. REFERENCES

1. Kamboj P.C., Pharmaceutical Analysis, Volume - I, Third edition, 2014, Vallabh Publications, Pages: 53–133.
2. Mendham J. et.al. Vogel's Textbook of Quantitative Chemical Analysis, Sixth edition, 2006, Dorling Kindersley (India) Pvt. Ltd, Pages: 33 - 36.

3. MATERIALS REQUIRED

Spatula, Funnel, Burette, Conical flask, Burette stand, Dropper, Volumetric pipette, Volumetric flask, Glass rod, Beaker, weighing balance.

4. THEORY

Terms involved in volumetric analysis:

1. **Titrant:** Solution whose concentration is known.
2. **Titrand:** Solution whose concentration is unknown.
3. **Stoichiometric / End point:** It shows that reaction between titrant and titrand is complete.
4. **Standard solution:** Solution whose exact concentration is known.
5. **Titration:** Reaction between titrand and titrant.

Steps involved in volumetric analysis:

1. **Method selection:** For analysis of base (NaOH) acid is used (HCl): Acid – Base titration.
2. **Sampling:** Small amount of chemical is taken as sample.
3. **Solution preparation:** Using appropriate formula, weight of chemical is calculated, weighed and dissolved in suitable solvent.
4. **Removing interferences:** Calibration, blank titration, parallel determination is done.
5. **Observation:** Volume of Titrant used for end point is observed.
6. **Calculation:** Using equivalent factor, concentration of sample solution is calculated.
7. **Result analysis:** Sample pass or fail as per pharmacopoeial standards

5. RESULTS

Steps involved in volumetric analysis were studied.

Signature of teacher:.....

EXPERIMENT NO. 2

1. OBJECTIVE

To prepare and standardize 100 ml of 0.1 N NaOH solution using Oxalic Acid as primary standard.

2. REFERENCES

1. Mendham J. et.al. Vogel's Textbook of Quantitative Chemical Analysis, Sixth edition, 2006, Dorling Kindersley (India) Pvt. Ltd, Pages: 40 – 59.

2. Kamboj P.C., Pharmaceutical Analysis, Volume - I, Third edition, 2014, Vallabh Publications, Pages: 270 – 271.

3. MATERIALS REQUIRED

Apparatus Required:

Funnel, volumetric flasks, conical flask, volumetric pipette, beaker, burette, burette stand, weighing balance.

Chemicals Required:

NaOH, oxalic acid, phenolphthalein solution.

4. THEORY

Primary standard is a chemical which has following properties: stable, high purity, non hygroscopic, high molecular weight and give stoichiometric reaction. Eg. oxalic acid, KCl.

Secondary standard is a chemical which has following properties: unstable, low purity, hygroscopic, low molecular weight, poor solubility. Eg. HCl, NaOH.

	Property	Primary Standard	Secondary Standard
1.	Purity	High	Low
2.	Stability	Stable	Unstable
3.	Oxidisation by air	No	Yes
4.	Hygroscopic	No	Yes
5.	Molecular weight	High	Low
6.	Solubility	Good	Poor

5. PROCEDURE

1. Preparation of NaOH solution

..... g of NaOH was weighed, dissolved in distilled water & volume made upto 100 ml.

2. Preparation of 0.1 N Oxalic acid solution.

.....g of Oxalic Acid was weighed, dissolved in distilled water & volume made upto 100 ml.

3. Standardization of NaOH solution

10 ml of Oxalic Acid solution & 1-2 drops of phenolphthalein solution were taken in a conical flask. To that solution NaOH solution was added until the solution became just pink.

The titration was performed 3 times along with 3 blank titrations. Volume of NaOH solution used in blank titration was subtracted from the titration value of NaOH to obtain the volume of NaOH solution used against Oxalic acid.

6. OBSERVATIONS

Table 1: Titration of 0.1 N Oxalic Acid solution Vs NaOH solution

S. No	Volume of 0.1 N Oxalic Acid solution (ml)	Burette Reading for NaOH Solution (ml)		
		Initial Reading	Final Reading	Used volume of NaOH
1.	10			(a)
2.	10			(b)
3.	10			(c)
Average volume of NaOH solution used in titration = $\frac{a+b+c}{3}$ = A				

Table 2: Blank Titration for 0.1 N Oxalic Acid solution Vs NaOH solution

S. No	Volume of Distilled water (ml)	Burette Reading for NaOH Solution (ml)		
		Initial Reading	Final Reading	Used volume of NaOH
1.	10			(x)
2.	10			(y)
3.	10			(z)
Average volume of NaOH solution used in blank titration = $\frac{x+y+z}{3}$ = B				

7. CALCULATIONS

Step 1: Calculations for weighing

$$\text{weight of the solute (w)} = \frac{N \cdot EW \cdot V}{1000}$$

1: 0.1 N NaOH solution

N = Normality of NaOH solution = 0.1

EW = Equivalent weight of NaOH = Molecular weight /1 =

V = Volume of NaOH solution = 100 ml

Thus w of NaOH taken = * * / 1000 = g

2: 0.1 N Oxalic Acid solution

N = Normality of Oxalic Acid solution = 0.1

EW = Equivalent weight of Oxalic Acid = Molecular weight /2 =

V = Volume of Oxalic Acid solution = 100 ml

Thus w of Oxalic Acid taken = * * / 1000 = g

Step 2: Calculation for Standardization of NaOH solution

$$N_{\text{NaOH}} V_{\text{NaOH}} = N_{\text{Oxalic Acid}} V_{\text{Oxalic Acid}}$$

Here

V_{NaOH} = volume used in titration (A) - volume used in blank titration (B) =

$N_{\text{Oxalic Acid}} = 0.1$

$V_{\text{Oxalic Acid}} = 10$

Thus
$$N_{\text{NaOH}} = \frac{N_{\text{Oxalic Acid}} V_{\text{Oxalic Acid}}}{V_{\text{NaOH}}} = \frac{0.1 \times 10}{\dots\dots\dots} = \dots\dots\dots \text{ N}$$

8. RESULTS

NaOH solution was prepared & was found to be N.

Signature of teacher:.....

EXPERIMENT NO. 3

1. OBJECTIVE

To perform the assay of ammonium chloride as per Indian Pharmacopoeia 2014.

2. REFERENCE

Indian Pharmacopoeia, Volume II, 2014, Pages: 1049-1050.

3. MATERIALS REQUIRED

Apparatus Required:

Funnel, volumetric flasks, conical flask, volumetric pipette, graduated pipette, burette, burette stand, rubber bulb, weighing balance.

Chemicals Required:

Sample (Ammonium chloride), formaldehyde solution, phenolphthalein solution, 0.1 M sodium hydroxide solution.

4. THEORY

“Acid base indicator is a chemical which determine the end point in neutralization titration by changing colour according to pH of titrand solution”.


Ostwald theory: Acid base indicator is either a weak acid or a weak base having different colour in ionized and unionized forms.

Example: Phenolphthalein (HIn) is a weak acid. It show pink coloured in basic medium and is colourless in acidic medium.



Phenolphthalein in basic solution: OH^- combine with H^+ of indicator to form H_2O . To maintain equilibrium, HIn is ionized (forward reaction increases) to increase the amount of In^- and indicator is colourless.

Phenolphthalein in acidic solution: HCl decreases the ionization of HIn due to common ion effect. This increases the rate of backward reaction (amount of unionized HIn increases) and indicator show pink colour.



5. PROCEDURE (I.P. 2014):

Dissolved 0.1 g ammonium chloride in 20 ml of water and added a mixture of 5 ml of formaldehyde solution, previously neutralized to dilute phenolphthalein solution, and 20 ml of water. After 2 minutes, titrated slowly with 0.1 M sodium hydroxide using a further 0.2 ml of phenolphthalein solution as indicator.

6. OBSERVATIONS

Table 1: Titration of ammonium chloride Vs NaOH solution

S. No	ammonium chloride (g)	Burette Reading for NaOH Solution (ml)		
		Initial Reading	Final Reading	Used volume of NaOH
1.	10			(a)
2.	10			(b)
3.	10			(c)
Average volume of NaOH solution = $\frac{a+b+c}{3}$ =				A

Table 2: Blank titration for ammonium chloride Vs NaOH solution

. No	Ammonium chloride (g)	Burette Reading for NaOH Solution (ml)		
		Initial Reading	Final Reading	Used volume of NaOH
1.	0			(x)
2.	0			(y)
3.	0			(z)
Average volume of NaOH solution = $\frac{x+y+z}{3}$ =				B

7. CALCULATIONS

Estimation of ammonium chloride

Equivalent factor: 1 ml of 0.1 M NaOH is equivalent to 0.005349 g of NH₄Cl.

Here

V_{NaOH} = volume used in titration (A) - volume used in blank titration (B) =

M_{NaOH} =

Each ml of 0.1 M NaOH is equivalent to

0.005349 g of NH₄Cl.

Thus each ml of M_{NaOH} NaOH is equivalent to

$M_{\text{NaOH}} \times 0.005349$ g of NH₄Cl.

Thus V_{NaOH} ml of M_{NaOH} NaOH is equivalent to

$V_{\text{NaOH}} \times M_{\text{NaOH}} \times 0.005349$ g of NH₄Cl.

= g of NH₄Cl.

= Experimental value

Claimed amount of NH_4Cl = quantity of NH_4Cl taken in the experiment = 0.1 g

Percentage Purity of ammonium chloride = $\frac{\text{Experimental amount of } \text{NH}_4\text{Cl}}{\text{Claimed Amount of } \text{NH}_4\text{Cl}} \times 100$

Claimed Amount of NH_4Cl

= _____ X 100 = %

8. RESULTS

The percentage purity of the given ammonium chloride was found to be %.

Signature of teacher:.....

EXPERIMENT NO. 4

1. OBJECTIVE

To prepare and standardize 100 ml of 0.05 N sodium thiosulphate solution using Potassium Iodate as primary standard.

2. REFERENCES

1.. Mendham J. et.al. Vogel's Textbook of Quantitative Chemical Analysis, Sixth edition, 2006, Dorling Kindersley (India) Pvt. Ltd, Pages: 456-466.

2.. Beckett A.H., Stenlake J.B., Practical Pharmaceutical Chemistry, Fourth edition - Part One, CBS Publishers. Page: 187.

3. MATERIALS REQUIRED

Apparatus Required:

Funnel, volumetric flasks, iodine flask, volumetric pipette, graduated pipette, burette, burette stand, rubber bulb, weighing balance.

Chemicals Required:

$\text{Na}_2\text{S}_2\text{O}_3$, KIO_3 , KI , H_2SO_4 .

4. THEORY

“**Iodometry** is a redox titration in which iodine is formed by reaction of sample with excess of iodide”. It is used in for quantitative reduction of oxidising agents. It involves:



In presence of light I^- is oxidised by oxygen to I_2 . Thus iodometry is done in diffused light / dark.



Titrand / titrant: Iodide (I^-). **Molecular iodine** is formed. Molecular iodine is highly volatile and sparingly soluble in water and it is retained in the titrand mixture due to formation of potassium tri-iodide with KI . **Indicator:** Iodine itself act as indicator (self indicator). Iodine is pale yellow in colour and iodide is colourless. Starch act as external indicator and iodine forms blue coloured complex with starch. **Iodine flask:** Iodometry is carried in iodine flask to prevent the Oxidation of KI by atmospheric oxidation & Loss of iodine formed due to the chemical reaction. **Oxidizing agent:** Example: KIO_3 , **Reducing agent:** I^- . **Titration:** I^- reacts with KIO_3 to form I_2 .

5. PROCEDURE

1. Preparation of $\text{Na}_2\text{S}_2\text{O}_3$ solution

.....g of $\text{Na}_2\text{S}_2\text{O}_3$ was weighed, dissolved in distilled water & volume made upto 100 ml.

2. Standardization of $\text{Na}_2\text{S}_2\text{O}_3$ solution

a. Preparation of 0.05 N KIO₃ solution.

...g of KIO₃ was accurately weighed, dissolved in distilled water & volume made upto 100 ml.

b. Titration of 0.05 N KIO₃ solution Vs Na₂S₂O₃ solution

10 ml of KIO₃ solution, 5.2 g of KI & 4 ml of 1 M H₂SO₄ were taken in an Iodine flask. Flask was closed, shaken & placed in dark for 10 minutes. To that solution Na₂S₂O₃ solution was added until the solution became colourless.

Blank titration can not be performed as in absence of KIO₃, iodine will not evolve to give colour to the titrand solution.

6. OBSERVATIONS**Table 1: Titration of Na₂S₂O₃ solution Vs 0.05 N KIO₃ solution**

S. No.	Volume of 0.05 N KIO ₃ solution (ml)	Burette Reading for Na ₂ S ₂ O ₃ solution (ml)		
		Initial Reading	Final Reading	Used volume of Na ₂ S ₂ O ₃
1.	10			(a)
2.	10			(b)
3.	10			(c)
Average volume of Na ₂ S ₂ O ₃ solution = $\frac{a+b+c}{3}$ =				

7. CALCULATIONS**Step 1: Calculations for weighing**

$$\text{weight of solute (g)} = \frac{N \cdot EW \cdot V}{1000}$$

1: 0.05 N Na₂S₂O₃ solution

$$N = \text{Normality of Na}_2\text{S}_2\text{O}_3 \text{ solution} = 0.05$$

$$EW = \text{Equivalent weight of Na}_2\text{S}_2\text{O}_3 = 158 / 1 = 158$$

$$V = \text{Volume of Na}_2\text{S}_2\text{O}_3 \text{ solution in ml} = 100 \text{ ml}$$

$$\text{Thus } w \text{ of Na}_2\text{S}_2\text{O}_3 \text{ taken} = * * / 1000 = \dots\dots\dots \text{ g}$$

2: 0.05 N KIO₃ solution

$$N = \text{Normality of KIO}_3 \text{ solution} = 0.05$$

$$EW = \text{Equivalent weight of KIO}_3 = 214 / 6 = 35.667$$

$$V = \text{Volume of KIO}_3 \text{ solution in ml} = 100 \text{ ml}$$

$$\text{Thus } w \text{ of KIO}_3 \text{ taken} = * * / 1000 = \dots\dots\dots \text{ g}$$

Step 2: Calculation for Standardization of Na₂S₂O₃ solution

$$N_{\text{Na}_2\text{S}_2\text{O}_3} V_{\text{Na}_2\text{S}_2\text{O}_3} = N_{\text{KIO}_3} V_{\text{KIO}_3}$$

Here $N_{\text{KIO}_3} = 0.05$

$$V_{\text{KIO}_3} = 10$$

$$V_{\text{Na}_2\text{S}_2\text{O}_3} =$$

$$\text{Thus } N_{\text{Na}_2\text{S}_2\text{O}_3} = \frac{N_{\text{KIO}_3} \times V_{\text{KIO}_3}}{V_{\text{Na}_2\text{S}_2\text{O}_3}}$$

=

8. RESULTS

Sodium thiosulphate solution was prepared & found to be ____N.

Signature of teacher:.....

EXPERIMENT NO. 5

1. OBJECTIVE

To prepare and standardize 100 ml of 0.02 N KMnO_4 solution using sodium oxalate as primary standard.

2. REFERENCES

1. Mendham J. et.al. Vogel's Textbook of Quantitative Chemical Analysis, Sixth edition, 2006, Dorling Kindersley (India) Pvt. Ltd, Pages: 447-450.

2. Beckett A.H., Stenlake J.B., Practical Pharmaceutical Chemistry, Fourth edition - Part One CBS Publishers. Pages: 177 – 178.

3. MATERIALS REQUIRED

Apparatus Required:

Funnel, volumetric flasks, conical flask, volumetric pipette, graduated pipette, burette, burette stand, rubber bulb, water bath, tripod stand, weighing balance.

Chemicals Required:

KMnO_4 , $\text{Na}_2\text{C}_2\text{O}_4$, sulphuric acid.

4. THEORY

Redox titration is based on change in electrode potential of sample during redox reaction. Oxidation is the loss of electrons / increase in oxidation number. Reduction is the gain of electrons / decrease in oxidation number.

Electrode potential is the tendency of element for loss or gain of electrons, when element is in contact with its ion. It is calculated by Nerst equation


$$E = E^\circ + \frac{0.059}{n} \log \frac{[\text{oxidized form}]}{[\text{reduced form}]}$$

E = electrode potential

E° = standard electrode potential

n = number of electrons gained or lost.

Redox curve: is curve between electrode potential of titrand and volume of titrant. Example: Titration of iodine (titrand) and sodium thiosulphate (titrant).



5. PROCEDURE

1. Preparation of KMnO_4 solution

..... g of KMnO_4 was weighed, dissolved in distilled water & volume made upto 100 ml.

2. Preparation of 0.02 N $\text{Na}_2\text{C}_2\text{O}_4$ solution

..... g of $\text{Na}_2\text{C}_2\text{O}_4$ was accurately weighed, dissolved in distilled water & volume made upto 100 ml.

3. Standardization of KMnO_4 solution

10 ml of $\text{Na}_2\text{C}_2\text{O}_4$ solution, 2.5 ml of 1 M sulphuric acid were warmed upto $60-70^\circ\text{C}$ in a conical flask using water bath. To that solution, KMnO_4 solution was added until the solution became just pink.

The titration was performed 3 times along with 3 blank titrations. Volume of the KMnO_4 solution used in blank titration was subtracted from the titration value of KMnO_4 to obtain the volume of KMnO_4 solution used against H_2O_2 solution.

6. OBSERVATIONS

Table 1: Titration of 0.02 N $\text{Na}_2\text{C}_2\text{O}_4$ solution Vs KMnO_4 solution

. S. No	Volume of 0.02 N $\text{Na}_2\text{C}_2\text{O}_4$ solution (ml)	Burette Reading for KMnO_4 Solution (ml)		
		Initial Reading	Final Reading	Used volume of KMnO_4
1.	10			(a)
2.	10			(b)
3.	10			(c)
Average volume of KMnO_4 solution = $\frac{a+b+c}{3} = \mathbf{A}$				

Table 2: Blank titration for 0.02 N $\text{Na}_2\text{C}_2\text{O}_4$ solution Vs KMnO_4 solution

. S. No	Volume of Distilled water (ml)	Burette Reading for KMnO_4 Solution (ml)		
		Initial Reading	Final Reading	Used volume of KMnO_4
1.	10			(x)
2.	10			(y)
3.	10			(z)
Average volume of KMnO_4 solution = $\frac{x+y+z}{3} = \mathbf{B}$				

7. CALCULATIONS

Step 1: Calculations for weighing

weight of the solute (w) = $N \cdot EW \cdot V / 1000$

1: 0.02 N $\text{Na}_2\text{C}_2\text{O}_4$ solution

N = Normality of $\text{Na}_2\text{C}_2\text{O}_4$ solution = 0.02

EW = Equivalent weight of $\text{Na}_2\text{C}_2\text{O}_4$ = Molecular weight / 2 =

V = Volume of $\text{Na}_2\text{C}_2\text{O}_4$ solution = 100 ml

Thus w of $\text{Na}_2\text{C}_2\text{O}_4$ taken = * * / 1000 = g

2: 0.02 N KMnO_4 solution

N = Normality of KMnO_4 solution = 0.02

EW = Equivalent weight of KMnO_4 in acidic medium = Molecular weight / 5 =

V = Volume of KMnO_4 solution = 100 ml

Thus w of KMnO_4 taken = * * / 1000 = g

Step 2. Calculation for Standardization of KMnO_4 solution

$$N_{\text{KMnO}_4} V_{\text{KMnO}_4} = N_{\text{Na}_2\text{C}_2\text{O}_4} V_{\text{Na}_2\text{C}_2\text{O}_4}$$

Here V_{KMnO_4} = volume used in titration (A) - volume used in blank titration (B)

$$N_{\text{Na}_2\text{C}_2\text{O}_4} = 0.02$$

$$V_{\text{Na}_2\text{C}_2\text{O}_4} = 10$$

$$\text{Thus } N_{\text{KMnO}_4} = \frac{N_{\text{Na}_2\text{C}_2\text{O}_4} V_{\text{Na}_2\text{C}_2\text{O}_4}}{V_{\text{KMnO}_4}}$$

$$= \dots\dots N$$

8. RESULTS

KMnO_4 solution was prepared & found to beN.

EXPERIMENT NO. 6

1. OBJECTIVE

To determine the percentage purity of given H₂O₂ sample (20 Volume).

2. REFERENCES

1. Kushwaha S.P. and Gupta S.K., Pharmaceutical Analysis - I, First edition, 2018, R. Narain Publishers & Distributors, Pages: 140 - 141.
2. Indian Pharmacopoeia 2014, Volume II, Page: 1912.

3. MATERIALS REQUIRED

Apparatus Required:


Funnel, volumetric flasks, conical flask, volumetric pipette, graduated pipette, burette, burette stand, tile, rubber bulb, water bath, tripod stand, weighing balance.

Chemicals Required:

KMnO₄, H₂O₂ (20 volume), sulphuric acid.

4. THEORY

“Redox indicator is a chemical which indicates the end point in redox titration by virtue of different colours in oxidized and reduced forms”.



When $n = 1$ and

- i. Electrode potential of titrand is more than ($E^{\circ}_{\text{In}} + 0.06$) colour of indicator is of oxidized form
- ii. Electrode potential of titrand is less than ($E^{\circ}_{\text{In}} - 0.06$) colour of indicator is of reduced form.
- iii. Electrode potential of titrand is between ($E^{\circ}_{\text{In}} \pm 0.06$) colour of indicator is due to mixture of oxidized and reduced form. This range ($E^{\circ}_{\text{In}} \pm 0.06$) is known as potential range and during this range, colour of indicator is not clear.

Example: Ferroin: $E^{\circ}_{\text{In}} = 1.14$ V

Types of redox indicators:

1. Self indicator: Titrand or titrant act as indicator. Eg: KMnO₄, I₂

KMnO ₄	Mn ⁺⁺
pink	colourless
2. Internal indicator: is added in the titration vessel (conical flask). eg: Diphenylamine.
3. External indicator: is used outside the titration vessel (conical flask). Titrand is taken out of the conical flask and added to indicator. Eg: Starch (iodine gives blue colour with starch).

5. PROCEDURE (I.P. 1996)

1 ml of H₂O₂ (20 Volume) solution, 20 ml of 1 M sulphuric acid were taken in a conical flask. To that solution 0.02 M KMnO₄ solution was added until the solution became pink.

Titration was performed 3 times along with 3 blank titrations. Volume of the KMnO₄ solution used in blank titration was subtracted from the titration value of KMnO₄ to obtain the volume of KMnO₄ solution used against H₂O₂ solution.

6. OBSERVATIONS

Table 1: Titration of H₂O₂ solution Vs KMnO₄ solution

S. No	Volume of H ₂ O ₂ solution (ml)	Burette Reading for KMnO ₄ Solution (ml)		
		Initial Reading	Final Reading	Used volume of KMnO ₄
1.	01			(a)
2.	01			(b)
3.	01			(c)
Average volume of KMnO ₄ solution = $\frac{a+b+c}{3}$ =				A

Table 2: Blank titration for H₂O₂ solution Vs KMnO₄ solution

S. No	Volume of H ₂ O ₂ solution (ml)	Burette Reading for KMnO ₄ Solution (ml)		
		Initial Reading	Final Reading	Used volume of KMnO ₄
1.	0			(x)
2.	0			(y)
3.	0			(z)
Average volume of KMnO ₄ solution = $\frac{x+y+z}{3}$ =				B

7. CALCULATIONS

Estimation of H₂O₂ strength

Here V_{KMnO_4} = volume used in titration (A) - volume used in blank titration (B)

=

$$N_{\text{KMnO}_4} = \dots\dots M_{\text{KMnO}_4}$$

Each ml of 0.02 M KMnO₄ is equivalent to 0.001701 g of H₂O₂

Thus each ml of M_{KMnO_4} KMnO₄ is equivalent to $M_{\text{KMnO}_4} \times 0.001701$ g of H₂O₂

Thus V_{KMnO_4} ml of M_{KMnO_4} KMnO₄ is equivalent to $V_{\text{KMnO}_4} \times M_{\text{KMnO}_4} \times 0.001701$ g of H₂O₂

= g of H₂O₂

= Experimental value

Claimed amount of H₂O₂ (quantity of H₂O₂ taken in the experiment) =

Percentage Purity of ammonium chloride = $\frac{\text{Experimental amount of H}_2\text{O}_2}{\text{Claimed Amount of H}_2\text{O}_2} \times 100$

= _____ X 100 = %

8. RESULTS

The percentage purity of the given H₂O₂ sample was found to be %.

EXPERIMENT NO. 7

1. OBJECTIVE

To prepare and standardize 100 ml 0.05 M disodium edetate solution.

2. REFERENCE

1. Kushwaha S.P. and Gupta S.K., Pharmaceutical Analysis - I, First edition, 2018, R. Narain Publishers & Distributors, Pages: 118 - 119.

2. Beckett A.H., Stenlake J.B., Practical Pharmaceutical Chemistry, Fourth edition - Part One CBS Publishers. Pages: 219 – 220.

3. MATERIALS REQUIRED

Apparatus Required:

Funnel, volumetric flasks, conical flask, volumetric pipette, graduated pipette, burette, burette stand, tile, rubber bulb, water bath, tripod stand, weighing balance.

Chemicals Required:

Disodium edetate, Mg ribbon, HCl, ammonia buffer solution, Eriochrome Black T, Sodium Chloride.


4. THEORY

“Complexometric titration is an analytical method in which volume of complexing agent used to convert sample into complex is measured”.

“Complexing agent (ligand) is an ion or molecule which form covalent or coordinate bond with metal ion to form complex”. It has O, N, S and donate electrons to metal ions to form coordinate bond.

Buffer solution

In complexometric reaction, acid is formed and reaction is reversible.



Complex of alkaline metal are stable in alkaline pH and decompose in neutral and basic medium.

Alkaline buffer are used to maintain pH 10. Buffer also prevent the backward reaction by neutralizing the acid formed in the titration reaction. Eg: NH_3 / NH_4Cl buffer.

4. PROCEDURE

1. Preparation of disodium edetate solution

..... g of disodium edetate was weighed, dissolved in distilled water & volume made upto 100 ml.

2. Standardization of disodium edetate solution

Weighted exactly about 0.24 g of metallic mg ribbon. dissolved metal in 25 ml of 1M HCl solution and diluted in volumetric flask to 1000 ml.

Transferred 25 ml of magnesium solution to conical flask, added 75 ml of distilled water, 2 ml of pH10 ammonia buffer solution and a pinch of eriochrome black T [1% (w/w) in sodium chloride].

Titrate with EDTA solution till the color changes to blue.

The titration was performed 3 times along with 3 blank titrations. Volume of the disodium edetate solution used in blank titration was subtracted from the titration value of disodium edetate to obtain the volume of disodium edetate solution used against Mg.

6. OBSERVATIONS

Step 1: Observation Tables

Table 1: Titration of Zn solution Vs disodium edetate solution

S. No	Volume of Zn solution (ml)	Burette Reading for disodium edetate Solution (ml)		
		Initial Reading	Final Reading	Used volume of EDTA
1.	10			(a)
2.	10			(b)
3.	10			(c)
Average volume of disodium edetate solution = $\frac{a+b+c}{3}$ =				A

Table 2: Blank titration for Zn solution Vs disodium edetate solution

S. No	Volume of Distilled water (ml)	Burette Reading for disodium edetate Solution (ml)		
		Initial Reading	Final Reading	Used volume of EDTA
1.	10			(x)
2.	10			(y)
3.	10			(z)
Average volume of disodium edetate solution = $\frac{x+y+z}{3}$ =				B

7. CALCULATIONS

Step 1: Calculations for weighing

weight of the solute (w) = $M \cdot MW \cdot V / 1000$

1: 0.05 M disodium edetate solution

M = Molarity of disodium edetate solution = 0.05

MW = Molecular weight of disodium edetate =

V = Volume of disodium edetate solution = 100 ml

Thus w of disodium edetate taken = * * /1000 = g

Step 3: Calculation for Standardization of disodium edetate solution

Equivalent factor: 1 ml of 0.1 M EDTA 0.638 mg of Zn

Here

$V_{\text{disodium edetate}} = \text{volume used in titration (A)} - \text{volume used in blank titration (B)} =$

Weight of Zinc in 10 ml = 0.024 g

Molarity of EDTA = weight of Zn in 10 ml X 1000

65.38 X volume of EDTA solution used in titration
= M

8. RESULTS

Disodium edetate solution was prepared & found to beM.

EXPERIMENT NO. 8

1. OBJECT

To determine the normality of base by conductometric titration.

2. REFERENCES

1. Kushwaha S.P. and Gupta S.K., Pharmaceutical Analysis - I, First edition, 2018, R. Narain Publishers & Distributors, Pages: 156.

2. Nema R.K., et. al., A Practical approach to Pharmaceutical Analysis, First edition, 2008, CBS Publisher, Pages: 116 - 123.

3. MATERIALS REQUIRED

Apparatus Required:

Funnel, beaker, glass rod, burette, burette stand, rubber bulb, conductometer.


Chemicals Required:

HCl sample solution, standard NaOH solution

4. THEORY

Conductance is the ability of ions to conduct electricity in an electrolyte solution. It is the reciprocal of resistance. Unit of conductance = $1/\text{ohm} = \text{ohm}^{-1} = \text{mho} = \text{Siemens}$

During titration, an ion is substituted by other ion. Conductance is increased or decreased depending on the conductance of latter ion (more or less mobility than earlier ion). Titrant should be 10-20 times more concentrated than titrand to prevent the change in titrand volume.



Conductivity cell contains sample and represents an arm of the Wheatstone bridge. Electrodes are made up of platinum black coated platinum plates (area of 1 cm^2) and are kept at a distance of 1 cm from each other. Electrodes are connected through platinum wires to Wheatstone bridge.

5. PROCEDURE

1. Filled the standard HCl solution in the burette
2. Took 25 ml of sample solution in a beaker.
3. Dipped the electrode in the sample solution.
4. After adjusting the cell constant, noted the conductance of the sample solution.
5. Added 1 ml of the standard HCl solution from the burette and stirred the sample solution for proper neutralization.
6. Noted the conductance of the sample solution again.

7. Repeatedly added 1 ml of standard HCl solution from the burette and stirred the sample solution to obtain readings.
8. Plotted the graph between conductance (y axis) and volume of standard HCl solution (x axis) added from the burette.


6. OBSERVATIONS

Step 1: Observation Table

Table 1: Titration of base vs acid by conductometric titration

S. No.	Volume of standard HCl solution added from the burette (ml)	Total Volume of standard HCl solution added from the burette (ml)	Conductance of sample solution
1.	0	0	
2.	1	1	
3.	1	2	
4.	1	3	
5.	1	4	
6.	1	5	
7.	1	6	
8.	1	7	
9.	1	8	
10.	1	9	
11.	1	10	

Step 2: Graphical Representation



7. CALCULATIONS

Normality of strong acid solution

$$N_{\text{HCl}} V_{\text{HCl}} = N_{\text{NaOH}} V_{\text{NaOH}}$$

Here

$$V_{\text{HCl}} = \text{volume used in titration} = \dots \text{ ml}$$

$$N_{\text{HCl}} = 1 \text{ N (given)}$$

$$V_{\text{NaOH}} = 25 \text{ ml}$$

$$\text{Thus } N_{\text{NaOH}} = \frac{N_{\text{HCl}} V_{\text{HCl}}}{V_{\text{NaOH}}} = \dots \text{ N}$$

8. RESULTS

Normality of base by conductometric titration was found to be $\dots \text{ N}$

EXPERIMENT NO. 9

1. OBJECTIVE

To determine the normality of strong base by potentiometric titration.

2. REFERENCES

1. Kushwaha S.P. and Gupta S.K., Pharmaceutical Analysis - I, First edition, 2018, R. Narain Publishers & Distributors, Pages: 166 - 167.

2. Nema R.K., et. al., A Practical approach to Pharmaceutical Analysis, First edition, 2008, CBS Publisher, Pages: 76 - 81.

3. MATERIALS REQUIRED

Apparatus Required:

Funnel, beaker, glass rod, burette, burette stand, rubber bulb, pH meter


Chemicals Required:

NaOH sample solution, standard HCl solution

4. THEORY

Glass electrode is used for determination of pH. It allow transfer of H^+ through membrane.

Glass electrode has bulb made up of thin glass membrane having H^+ ion permeability and high electrical conductivity.



Glass electrode develops a potential at the membrane separating the standard HCl and the test solutions due to H^+ concentration gradient. Glass bulb function act as semipermeable membrane to H^+ . Ion exchange reaction takes place between strongly charged cation in glass lattice and proton from the solution.

5. PROCEDURE

1. Filled the standard HCl solution in the burette
2. Took 25 ml of sample solution in a beaker.
3. Dipped the electrode in the sample solution.
4. Noted the electrode potential of the sample solution.

5. Added 1 ml of the standard HCl solution from the burette and stirred the sample solution for proper neutralization.
6. Noted the electrode potential of the sample solution again.
7. Repeatedly added 1 ml of standard HCl solution from the burette and stirred the sample solution to obtain readings.
8. Plotted the graph between electrode potential (y axis) and volume of standard HCl solution (x axis) added from the burette.


6. OBSERVATIONS

Table 1: Titration of base vs acid by potentiometric titration

S. No.	Volume of standard HCl solution added from the burette (ml)	Total Volume of standard HCl solution added from the burette (ml)	eV
1.	0	0	
2.	1	1	
3.	1	2	
4.	1	3	
5.	1	4	
6.	1	5	
7.	1	6	
8.	1	7	
9.	1	8	
10.	1	9	
11.	1	10	

Step 2: Graphical Representation

Object: To determine the normality of base solution by potentiometric titration.



7. CALCULATIONS

Normality of strong base solution

$$N_{\text{NaOH}} V_{\text{NaOH}} = N_{\text{HCl}} V_{\text{HCl}}$$

Here

$$V_{\text{HCl}} = \text{volume used in titration} = \dots \text{ ml}$$

$$N_{\text{HCl}} = 1 \text{ N (given)}$$

$$V_{\text{NaOH}} = 25 \text{ ml}$$

$$\text{Thus } N_{\text{NaOH}} = \frac{N_{\text{HCl}} V_{\text{HCl}}}{V_{\text{NaOH}}} = \dots \text{ N}$$

8. RESULTS

Normality of base by potentiometric titration was found to be $\dots \text{ N}$

EXPERIMENT NO. 10

1. OBJECTIVE

To draw neutralization curve between the given acid-base samples & predict the indicator that can be used.

2. REFERENCES

1. Kushwaha S.P. and Gupta S.K., Pharmaceutical Analysis - I, First edition, 2018, R. Narain Publishers & Distributors, Pages: 166 - 167.

2. Nema R.K., et. al., A Practical approach to Pharmaceutical Analysis, First edition, 2008, CBS Publisher, Pages: 76 - 81.

3. MATERIALS REQUIRED

Apparatus Required:

Funnel, beaker, glass rod, burette, burette stand, rubber bulb, pH meter.

Chemicals Required:

NaOH solution, HCl solution

4. THEORY

“Acid base indicator is a chemical which determine the end point in neutralization titration by changing colour according to pH of titrand solution”.


Ostwald theory: Acid base indicator is either a weak acid or a weak base having different colour in ionized and unionized forms. Example: **Phenolphthalein** (HIn) is a weak acid. It show pink coloured in basic medium and is colourless in acidic medium.



Phenolphthalein in basic solution:

OH^- combine with H^+ of indicator to form H_2O . To maintain equilibrium, HIn is ionized (forward reaction increases) to increase the amount of In^- and indicator is colourless.

Phenolphthalein in acidic solution:



HCl decreases ionization of HIn due to common ion effect. This increases the rate of backward reaction (amount of unionized HIn increases) and indicator show pink colour.

5. PROCEDURE

1. Filled the NaOH solution in the burette
2. Took 25 ml of HCl solution in a beaker.
3. Dipped the electrode in the HCl solution.
4. Noted the pH of the HCl solution.


5. Added 1 ml of the NaOH solution from the burette and stirred the HCl solution for proper neutralization.
6. Noted the pH of the HCl solution again.
7. Repeatedly added 1 ml of NaOH solution from the burette and stirred the HCl solution to obtain readings.
8. Plotted the graph between pH (y axis) and volume of NaOH solution (x axis) added from the burette.

6. OBSERVATIONS

Table 1: Titration of base vs acid by pH meter

S. No.	Volume of NaOH solution added from the burette (ml)	Total Volume of NaOH solution added from the burette (ml)	pH
1.	0	0	
2.	1	1	
3.	1	2	
4.	1	3	
5.	1	4	
6.	1	5	
7.	1	6	
8.	1	7	
9.	1	8	
10.	1	9	
11.	1	10	

Step 2: Graphical Representation



7. RESULTS

The neutralization curve between the given acid-base samples was drawn & the indicator\ s that can be used is\ are.....

GLOSSARY

- **Sample:** Small amount of pharmaceutical substance representing the property of complete pharmaceutical substance.
- **Analyte:** The component of interest in the sample.
- **Impurity:** A component present in the pharmaceutical product that is not the chemical entity defined as the drug substance or an excipient in the drug product.
- **Pharmaceutical analysis:** The identification and chemical composition determination of the analyte present in the sample of the pharmaceutical substances.
- **Quantitative analysis:** It determine the amount of analyte and impurity present in the sample.
- **Molarity (M):** It indicates the number of moles of solute present per litre of solution". Its unit is moles/litre.
- **1 mole:** Total of the atomic weights of the elements present in the solute. Example: 1 mole NaOH = 23 +16 + 1 = 40g.
- **Molar solution (1M):** It contains 1 mole of solute present in 1 litre of solution. Example: 1 M NaOH solution contains 40 g NaOH in 1 litre of solution.
- **Normality (N):** It indicates the number of gram equivalents of solute present per litre of solution". Its unit is gm equivalents/litre.
- **Volumetric titration:** Reagent solution is added to sample solution (or vice versa) until reaction is complete. Volume of reagent used to achieve end point of the chemical reaction is measured and used for calculating the amount of analyte present in the sample.
- **Titration:** Reaction between titrand and titrant.
- **Standard solution:** It is a solution whose exact concentration is known. It is known as titrant during titration.
- **Titrand:** It is a solution used in titration whose exact concentration is unknown.
- **Standardisation:** It is the process of determining the exact concentration of a solution.
- **Stoichiometric point / End point / equivalence point:** It is the titration point at which chemical reaction between titrant and titrand is complete and titrant and titrand are in equilibrium. Eg: At end point in acid – base titration, moles of H⁺ are equal to moles of OH⁻
- **Indicator:** is a chemical which indicates the end point of the titration.
- **Acid base (Neutralization) titration:** Sample / reagent are complimenting acid and base and undergo neutralization reaction with each other.

-
- **Oxidation-reduction (Redox) titration:** Sample / reagent are complimenting oxidizing agent and reducing agent and undergo redox reaction with each other.

- **Complexometric Titration:** Sample and reagent react to form a complex.


$$\text{CaCO}_3 + \text{EDTA} \rightarrow \text{Ca-Edetate}$$

- **Electrochemical methods of analysis:** They determine the end point in a titration by measuring the change in electrical property of the sample during the chemical reaction.

- **Conductometry:** It measure the change in conductance due to substitution of one ion by other ion.

- **Potentiometry:** It measure change in electrode potential of sample during redox reaction.

- **Redox reaction:** It is a reaction in which oxidation and reduction takes place together”.



- **Oxidizing agent (oxidant):** It help in oxidation of other chemical and itself gets reduced. Eg. Ceric sulphate, I_2 , KIO_3 .

- **Reducing agent (reductant):** It help in the reduction of other chemical and itself gets oxidized. Eg: KI , $\text{Na}_2\text{S}_2\text{O}_3$.

- **Neutralization curve:** It is a curve between pH of titrand (x axis) and volume of titrant used in acid base titration (y axis).